

The only software that provides automated and objective assessment of exam quality on every mammogram.

Smart dashboards, alerts, and reports optimize breast cancer screening operations. Volpara® Analytics™ software makes it easier to provide evidence to the ACR and FDA that centers are meeting quality standards. This vendor-neutral software is used by facilities around the world to help ensure they provide high-quality care to the women in their communities.



THE HEART OF ANALYTICS

Volpara® Science™

Our clinically validated algorithms that use x-ray physics and artificial intelligence to assess image quality and breast tissue composition.

300+

peer-reviewed articles and research abstracts that validate Volpara technology



TruPGMI™

Automatically and objectively assess technologists' patient positioning and resulting image to help make sure that every mammogram is of diagnostic quality. It provides technologists with fair performance feedback.



TruRadDose™

Analyse the radiation dose delivered to patients based on their breast density instead of the dose estimated by the equipment manufacturer.



TruPressure™

Compute the compression pressure applied to each patient's breasts. This helps technologists better understand the patient experience and improves the clinical effectiveness of the mammogram.




TruDensity™

Automatically assess volumetric breast density on a continuous scale. Its precise, objective values help radiologists provide personalized screening and understand the sensitivity of each mammogram.

Streamlined compliance and accreditation prep

> FDA MQSA EQUIP

Use Volpara Analytics for easy FDA MQSA EQUIP compliance preparation. Rapidly find and fix image quality issues. Schedule EQUIP reports to generate with your objective Volpara quality data to show inspectors.



Address the required inspection questions:

①

Corrective actions

Use comprehensive objective trend reporting to identify areas for staff improvement.

②

Periodic image quality review

Generate detailed reports to show inspectors.

③

Quality control oversight

Demonstrate lead interpreting physician oversight of QA/QC and corrective action processes.

> ACR accreditation

Volpara Analytics simplifies your ACR accreditation efforts, reducing the time spent searching PACS for your technologists' best work. Lead technologists can simply set the date range and filter by machine, TruPGMI Perfect or Good positioning scores, and breast density category. Volpara Analytics then produces a short list of cases that's easy to review.

“I can't see how a large and busy practice in breast imaging can meaningfully meet quality assurance without Volpara Analytics.”


– Robin Shermis, MD, MPH, FACR, Medical Director, ProMedica Breast Care

Continuous quality improvement

> Fair performance review

Your technologists work hard and want to feel that their performance is assessed fairly. Volpara Analytics provides objective feedback to inform performance reviews.

Technologists can log in and immediately see how to improve their skills and image quality so they can proactively focus on improving their performance over time. They can trust that their feedback is not based on a handful of cases. Every single image is evaluated. Lead technologists, rather than being restricted to biannual assessments, can easily review performance with each technologist at any time.



40 million+

mammography images have been evaluated for positioning and compression by Volpara

3,600+


technologists use Volpara to monitor performance

2,000+

facilities have installed Volpara software, including top US cancer centers

> Empowered technologist training

Using Volpara Analytics as a personalized coach, technologists have on-demand access to an extensive library of positioning videos and guides produced in collaboration with Mammography Educators and noted expert Louise Miller, R.T. (M)(ARRT), CRT(M), FSBI, FNCBC.



> Building a quality culture

The quality metrics in Volpara Analytics create an environment of continuous improvement without placing undue burden on lead technologists or supervisors. Volpara Analytics uniquely displays technologists' metrics alongside the corresponding images so they can easily visualize what they need to do better.

With so much information readily available, your team may also experience an internal shift as their day-to-day work is more directly connected to wider quality initiatives in your organization and promoting the facility's reputation and growth.

“Volpara Analytics contains the amazing capability to analyze patient positioning. We can now get objective measurements and feedback on the quality of positioning by our technologists. We can then personalize training based on their needs and partner them with other staff who excel in their area of opportunity.”

– Shakira Sarquis-Kolber, Director of Women's Imaging, Christine E. Lynn Women's Health & Wellness Institute, Boca Raton Regional Hospital

Business success story: Mayfair Diagnostics

Mayfair Diagnostics is a busy, multi-modality imaging practice in Calgary, Alberta, Canada. With 44 technologists across 10 breast imaging locations, Mayfair takes 54,000 mammograms annually.

Here's a look at the improvements Mayfair saw within the first six months of using Volpara Analytics.*

*Measures were calculated using PGMI prior to the May 2018 adoption of the new Volpara TruPGMI method of evaluating the clinical image quality of mammograms. Actual percentage improvements may vary with the new standard.

20%

increase in target compression

60%

reduction in time searching for cases

37%

reduction in inadequate images

14%

reduction in technical recalls

MAYFAIR[®]
DIAGNOSTICS

“There is something comforting about having the quality assessed on every mammogram I take. It gives me constant and very specific reminders about what I can do to improve and provides me with a rewarding feeling when I'm doing a good job. I can't imagine doing mammography without Volpara Analytics.”

– Katie Kaminski, Mammography Coordinator, Mayfair Diagnostics

Additional benefits

Volpara Analytics provides the secure,* actionable intelligence needed by your facility, large or small, to operate at peak performance:

Dynamic, interactive dashboards

Designed specifically for the major roles at your breast imaging facility—from a breast imaging manager studying physician referral patterns to a risk counselor identifying patients who may benefit from essential screening.

Vendor-neutral analytics

To monitor equipment utilization and exam time for resource planning.

Proactive dose monitoring

To help ensure that women get the appropriate dose calculation for their breast density—regardless of manufacturer.

*Volpara Analytics maintains patient data security through Volpara's ISO 27001:2013 certification. In addition, the Volpara Breast Health Platform benefits from Microsoft Azure's platform compliance.

Breast Health Platform

Volpara Analytics is only one part of Volpara's software platform, which can easily be scaled as centers are added, volume increases, or patient care needs evolve. Additional integrated modules—including software for lung screening programs—can be quoted and configured at your request.



[Analytics](#)



[Scorecard](#)



[Patient hub](#)



[Risk](#)



[Live](#)

13.5 million⁺

women across 39 countries have had their breast composition assessed by Volpara Health

With the most scientifically validated technology platform in the breast health industry, Volpara Health is the preferred partner of leading clinical sites around the world. Every day, Volpara remains focused on its mission to help save families from cancer.

Interested in Analytics?

Contact your Volpara Health representative for a live demonstration, or visit volparahealth.com.

Contact

info@volparahealth.com
support@volparahealth.com

US +1 855 607 0478
AUS 1800 370 623
NZ 0800 444 148

Europe +44 203 051 1029
Global +64 4 499 6029

Connect

 @VolparaHealth

 @volpara

 Volpara Health

References

1. Clinical outcomes in very early breast cancer ($\leq 1\text{cm}$): A national population based analysis. Mahvish Muzaffar, Abdul Rafeh Naqash, Nasreen A. Vohra, Darla K. Liles, and Jan H. Wong Journal of Clinical Oncology 2017 35:15-suppl, e12034-e12034

2. Holland, K., et al., Performance of breast cancer screening depends on mammographic compression, in Breast Imaging: 13th International Workshop, IWDM 2016, Malmö, Sweden, June 19–22, 2016, Proceedings, A. Tingberg, K. Lång, and P. Timberg, Editors. 2016, Springer International Publishing: Cham. p. 183–189.

3. Whelehan, P., et al., The effect of mammography pain on repeat participation in breast cancer screening: a systematic review. Breast, 2013. 22(4): p. 389–94.

Copyright © Volpara Health Technologies Ltd. 2021. All rights reserved.